

UCP1 - Mitochondrial brown fat uncoupling protein 1

Product specification

Acronym: UCP1	Purity: >60%
Class: Transporter	Activity: Proven
Origin: Mitochondrial	Length: Full Length
Molecular weight: 33 kDa	TMD: 6
Application: Screening & Display Technologies	Biological function: Thermogenesis

Product description

UCP1 is a mitochondrial transporter protein that create proton leaks across the inner mitochondrial membrane, thus uncoupling oxidative phosphorylation from ATP synthesis. As a result, energy is dissipated in the form of heat.

Protein Source: UCP1 wild type protein

Fig.1: AA sequence of UCP1 protein

10	20	30	40	50
MGGLTASDVH	PTLGVQLFSA	GIAACLADVI	TFPLDTAKVR	LQVQGECPST
60	70	80	90	100
SVIRYKGVLG	TITAVKTEG	RMKLYSGLPA	GLQRQISSAS	LRIGLYDTVQ
110	120	130	140	150
EFLTAGKETA	PSLGSKILAG	LTTGGVAVFI	GQPTEVVKVR	LQAQSHLHGI
160	170	180	190	200
KPRYTGTYNA	YRIIATTEGL	TGLWKGTTPN	LMRSVIINCT	ELVTYDLMKE
210	220	230	240	250
AFVKNILAD	DVPCHLVSAL	IAGFCATAMS	SPVDVVKTRF	INSPPGQYKS
260	270	280	290	300
VPNCAMKVFT	NEGPTAFFKG	LVPSFLRLGS	WNVIMFVCFE	QLKRELSKSR
QTMDCAT				

Affinity Tag: Histidine tag fused to the N-terminal end of the protein.


Production conditions: UCP1 is expressed in a cell-free expression system in the presence of lipid vesicles. 100 µg can be produced and qualified in about 1 week.

Quality analysis

Purity: Typically >60% as determined by SDS-Page and Coomassie Blue staining.

Purification procedure: As standard, UCP1 proteoliposomes are purified on a sucrose gradient. Further purification steps can be added if required.

Fig.2: Proteoliposome UCP1 after purification (Western blot identification).



Assessment of functionality

Cell-free expression systems provide a real alternative for membrane protein expression, enabling the study of structure and function of membrane proteins.

Methods: SPR ATP binding assay

Results:

The binding properties of UCP1 proteoliposomes have been validated using Horiba Scientific SPR ATP binding assay. Small molecule were injected on a biochip grafted with UCP1 biotinylated proteoliposomes. We detected specific interactions between UCP1 proteoliposomes and different concentrations of small molecules. The signal was dose dependent. Different liposome composition were tested. A positive and dose dependant signal was observed with the liposome composition D and the K_D was estimated to 350 μM .




Fig.3: SPR results. A /Sensorgrams obtained at different concentrations of small molecules and with different liposome composition. B/ Variations of reflectivity obtained at steady state. B) Reflectivity after negative control subtracting.

Formulation


Buffer: Available in Tris 50mM, pH 7.5. Other buffers or customized formulation can be provided upon request.

Customized Hydrophobic matrix: Customized formulation with specific lipids like PEGylated or biotinylated lipids can be used upon request, as well as targeting molecules.

Storage/Stability: Store at +4°C for up to one week or several months at -80°C. Aliquot for storage. Do not freeze-thaw after aliquoting.

Use restrictions: For life science research use only.

Available sizes: 10µg, 20µg, 100 µg, 200 µg, 500 µg, bulk



Need a specific amount, a quote or
any additional information?
Contact-us



Synthelisis

T : +33 (0)4 76 54 95 35
E: contact@synthelisis.fr
www.synthelisis.com

